

**Patterson® PDCARE™ Wipes and PDCARE™ Surface
(087-7191, 087-7209, 085-0602, 085-0610, and 085-0628)
Disinfectant Kill Times**

- 3 Minutes
 - *Mycobacterium tuberculosis* var: bovis (BCG)
 - *Staphylococcus aureus*
 - *Pseudomonas aeruginosa*
 - *Salmonella enterica*
 - *Trichophyton mentagrophytes*
- 2 minutes
 - Methicillin Resistant *Staphylococcus aureus* (MRSA)
 - Vancomycin Resistant *Enterococcus faecalis* (VRE)
 - *Staphylococcus aureus* with Reduced Susceptibility to Vancomycin
 - Hepatitis B Virus (HBV)
 - Hepatitis C Virus (HCV)
 - Herpes Simplex Virus Type 1
 - Herpes Simplex Virus Type 2
 - Human Immunodeficiency Virus (HIV-1)
 - Influenza A2 Virus