



G-aenial™ Universal Flo


UNIVERSAL INJECTABLE COMPOSITE RESTORATIVE



THE FIRST TRULY UNIVERSAL INJECTABLE COMPOSITE



Results of G-aenial Universal Flo Evaluation
at 1 Year from The Dental Advisor



TOUCH of GENIUS™

- Superior polish and high gloss retention¹
- Lack of post-operative sensitivity
- Extraordinary wear resistance
- Strength through filler innovation - for Classes I, II, III, IV and V
- Wears like a traditional composite, yet handles like a flowable



The First Truly Universal Injectable Composite



Truly Universal - Indicated for:



Class I



Class II



Class III



Class IV



Class V


True innovation in filler size and composition, in filler silanation and in manufacturing techniques, have given rise to a remarkable composite resin. G-ænial Universal Flo is the first truly universal Injectible Composite, a material that clinicians can place with ease for beautifully adapted, beautifully finished, fast, aesthetic and tough composite restorations.

Restored using shades A03, A4, A3 and AE.

Dr. S. Koide, Nagano, Japan




BEFORE




20.0 µm

G-ænial Universal Flo - Features 200-nm strontium glass fillers homogeneously dispersed for high flexural strength and wear resistance



20.0 µm

Filtek Supreme Ultra Flow⁺ - (Nano composite). Typically contains spherical nano clusters made up of agglomerated nano fillers of varying sizes



20.0 µm

Clearfil Majesty⁺ - (Micro Hybrid). Typically contains a range of glass and nano fillers of varying sizes

400222 **G-ænial™ Bond & G-ænial™ Universal Flo Kit**
 Contains: One G-ænial Bond Kit (5 mL), 1 Syringe A2 G-ænial Universal Flo and 20 Dispensing Tips.

G-ænial™ Universal Flo Refills
 Contains: One Syringe and 20 Dispensing Tips.
 004202 A1 004203 A2 004204 A3 004205 A3.5 004206 A4 004216 AE
 004207 B1 004208 B2 004209 B3
 004213 BW 004210 C3 004214 CV 004215 JE

G-ænial™ Flo Refills
 Contains: One Syringe and 20 Dispensing Tips.
 004194 A1 004195 A2 004196 A3 004197 A3.5 004198 A4 004199 AO2
 004200 AO3 004201 CV

004635 **Metal Dispensing Tips III (x30)**
 004634 **Plastic Dispensing Tips III (x30)**



¹Data on file.
[†]Not a registered trademark of GC America Inc.