

# Herculite® *Ultra Flow*

Nanohybrid Flowable Composite



A FLOWABLE THAT WORKS  
**JUST THE WAY YOU WANT.**

Your practice is *our* inspiration.™

# The Ultimate Flowable. The Ultimate Match.

## THE PERFECT DUO

For decades, the Herculite® brand has established clinical success in the practice. Now, Herculite Ultra Flow combines the long-standing expertise of the Herculite brand with a new innovative flowable composite. Plus, it's the perfect match with the esthetics and performance of Herculite Ultra.



Base/Liner

Case image courtesy of Dr. Bob Lowe.


## SMART PLACEMENT TECHNOLOGY

For a flowable, having one viscosity that flows into all internal line and point angles while holding its shape may seem impossible. However, the rheological additive in Herculite Ultra Flow helps it perform this feat by reducing surface tension as the material contacts the tooth surface under pressure – allowing it to flow to the boundaries of the restoration.

The same additive encourages higher surface tension after the material has been placed. So, it won't run or drip.


# A flowable with the performance you need.

## OUTSTANDING FLEXURAL STRENGTH\*




Similar to its companion Herculite Ultra, Herculite Ultra Flow has excellent mechanical properties, including strength.

## LOW SHRINKAGE\*




Leading experts cite low shrinkage as a critical feature of flowable composites; Herculite Ultra Flow was specially designed with a high-filler loading rate for minimal polymerization shrinkage.

## OVER 300% RADIOPACITY\*



With its high radiopacity, Herculite Ultra Flow is easily discernable on radiographs.

## HIGH GLOSS RETENTION\*



While initial gloss is a critical feature, gloss over time is more important. It is an indication of how long the esthetics of the restoration will last. Herculite Ultra Flow retains gloss after aging better than many other leading materials.

<sup>a</sup> Contrast gloss = Specular reflection @ 60° / Diffuse reflection

\*Data from Kerr R&D

Your practice is *our* inspiration.™

**Kerr**™

# ORDER TODAY!

For your convenience, Herculite Ultra Flow comes in 12 shades.

## REFILLS:

- 2 Syringes (2g each)
- 20 Dispensing Tips
- Technique Guide
- Directions for Use

## ACCESSORIES:

33370 - Tip Refills (pack of 50)



## SHADES:

- 35392 A1
- 35393 A2
- 35394 A3
- 35395 A3.5
- 35396 A4
- 35397 B1
- 35398 B2
- 35399 C2
- 35400 D2
- 35401 XL1
- 35402 XL2
- 35403 Universal Opaque

## RELATED PRODUCTS

### BONDING AGENTS

- OptiBond® Solo Plus™
- OptiBond® XTR
- OptiBond® FL
- OptiBond® All-In-One



### CURING LIGHTS

- Demi™ Ultra
- Demi™ Plus



### UNIVERSAL COMPOSITE

- Herculite® Ultra
- Premise™
- Herculite® XRV™
- Point 4™



### SINGLE FILL™

- SonicFill™



**Kerr™**