



Introducing the NEW E-STATIS™ family: the E-STATIS Advantage and E-STATIS Prep.




Products in the E-STATIS Family come equipped with the new and advanced E-STATIS SLM (Short, LED, Micromotor). The SLM is 35% shorter and lighter than today's current motors. The E-STATIS Advantage offers a range of customizable settings for both operative and endodontic procedures. The E-STATIS Prep offers customizable settings for up to 8 preparatory procedures. And last but not least, the remarkable touch screen control is detachable for modular placement.

Leave it to the infection control specialist to develop one of the most technologically advanced motors that is fully autoclavable!

For more information on the E-STATIS Advantage or E-STATIS Prep please visit www.scican.com.

Your Infection Control Specialist™

SciCan Dental



* only E-STAT/S Advantage

SPECIFICATIONS

ADVANTAGE w/ torque control

PREP

Motor speed range:	100 - 40,000 min-1	2,000 to 40,000 min-1
Power supply:	120 W	60 W
Motor voltage:	max. 22 V AC	max. 22 V AC
Torque:	max. 3.5 Ncm	max. 3 Ncm
Operation current of the LED:	max. 250 mA	max. 250 mA
Air outlet of coupling (cooling air)	3.5 - 8 NL/min	3.5 - 8 NL/min
Spray air pressure:	1 - 2.5 bar / 14.5 - 36.2 psi	1 - 2.5 bar / 14.5 - 36.2 psi
Spray water pressure:	0.8 - 2.0 bar / 11.6 - 29.0 psi	0.8 - 2.0 bar / 11.6 - 29.0 psi
Weight (motor only)	62g	62g
Dimensions control unit:	4.33 x 2.99 x 6.89 inch (110 x 76 x 175 mm)	4.33 x 2.99 x 6.89 inch (110 x 76 x 175 mm)
Dimensions motor:	1.25 x 0.87 inch (31.7 x 22 mm)	1.25 x 0.87 inch (31.7 x 22 mm)
Rotation:	clockwise / counter clockwise	
Connection of attachments:	ISO 3964	ISO 3964
# of Presets:	3 x prep, 5 x endo	8 x prep

PRODUCT CODE:

S21200	E-STAT/S Advantage Electric Handpiece System
S21300	E- E-STAT/S Prep Electric Handpiece System
S21103	E-STAT/S Display extension kit

Manufactured by:
SciCan GmbH
 Wangener Strasse 78
 88299 Leutkirch
 GERMANY
 Tel.: +49 (0)7561 98343 - 0
 Fax: +49 (0)7561 98343 - 699

SciCan Inc.
 701 Technology Drive
 Canonsburg, PA 15317
 USA
 Phone: 724 820 1600
 Toll free: 1-800 572 1211
 Fax: 724 820 1479

SciCan Ltd.
 1440 Don Mills Road
 Toronto ON M3B 3P9
 CANADA
 Phone: 416 445 1600
 Toll free: 1-800-667-7733
 Fax: 416 445 2727

Three year warranty.

Fully autoclavable.

Brushless technology – greatly reduces service requirements.

One of the shortest & lightest torque-controlled dental motors in the world – weighs only 62 grams.

Adjustable LED light for unsurpassed visibility



Compact design for easy installation.



Hose connector: ISO 9168 fits 4, 5 or 6 hole connection.



Extension kit provides flexibility to remotely locate both control box and touchpad.

