

SmartLite®  
**max**  
L.E.D. Curing Light



## SPACE-AGE TECHNOLOGY. NEW-AGE AFFORDABILITY.

### A ROBUST FEATURE SET. ALL IT NEEDS IS YOU.

The SmartLite® Max LED Curing Light provides an array of features so you can do your best work, conveniently and flexibly.

- **Dual Wavelength LEDs:** Be confident that any restorative material is being thoroughly cured the first time around. (Dual LEDs are for photoinitiators in the 377-490nm wavelength range.)
- **Powerful precision:** High output of up to 2850mW/cm<sup>2</sup> and a large curing area.
- **Multiple modes:** With boost, ramp, pulse and standard mode you can perform the way that suits each dental case best.
- **More mobility:** With a long-lasting battery, and cordless or corded operation, you can avoid running out of battery in the middle of a procedure.
- **Easy cleanup:** The SmartLite® Max LED Curing Light is fan- and vent-free to allow for easy disinfection.

# YOU CAN HAVE IT ALL.

It's not rocket science that better tools allow you to do your job more efficiently. And more affordability is just smart business. That's why we designed the SmartLite® Max LED Curing Light with dual wavelength technology, multiple modes, and mobility features — all at a down-to-earth price.



## DESIGN AND FEATURES THAT ARE LIGHT YEARS AHEAD.

The **SmartLite® Max Curing Light** is at the forefront of LED curing light technology. It features a sleek, hi-tech design and all the features you need to do your best work, conveniently and efficiently.

## SPACE-AGE TECHNOLOGY THAT OUTSHINES THE COMPETITION.

Brand	Sold By	Highest Claimed Output (mW/cm <sup>2</sup> )*	Multiple Wavelength	Use Corded or Cordless?	Ramp Mode	Boost Mode	Pulse Mode	Fan	Built-in Radiometer
<b>SmartLite® Max</b>	<b>DENTSPLY</b>	<b>2850</b>	✓	✓	✓	✓	✓	—	✓
Bluephase 20i	Vivadent	2000	✓	✓	✓	✓	—	✓	—
Coltolux LED	Coltene	1000	—	—	—	—	—	—	—
Demi	Kerr	1330	—	—	—	—	✓	✓	—
Elipar S10	3M	1200	—	—	—	—	—	—	✓
FLASHlite 1401	Discus	1100	—	—	—	—	—	—	✓
Radii-cal	SDI	1200	—	—	✓	—	✓	—	✓
Ultra-Lume 5	Ultradent	800	✓	—	—	—	—	—	—

Visit [www.smartlitemax.com](http://www.smartlitemax.com) or call your DENTSPLY Caulk rep for more information.

## ORDERING INFO

- 644050 **SmartLite Max Complete Unit 120v**  
1 Charging Base Assembly, 1 Handpiece Assembly,  
1 Curing Probe, 1 Wall Plug-in Power Supply 120v,  
1 Eye Protection Assembly, 50 Disposable Barrier Sleeves
- 644052 **SmartLite Max Power Supply 120v**  
1 Wall Plug-in Power Supply, Standard Adaptor
- 644054 **SmartLite Max Charging Base Assembly (1)**
- 644056 **SmartLite Max Replacement Curing Probe (1)**
- 644057 **SmartLite Max Eye Protection Assembly**  
1 Eye Protection Shield, 1 Grommet
- 644060 **SmartLite Max Disposable Barrier Sleeves (250)**

**DENTSPLY**  
CAULK

38 West Clarke Avenue  
Milford, DE 19963  
[www.dentsply.com](http://www.dentsply.com) • [www.caulk.com](http://www.caulk.com)  
1-800-LD-CAULK (532-2855)  
In Canada, call 1-800-263-1437