



**Product Name:** Patterson® pdCARE™ Dual Enzymatic Detergent, 1 Gallon

**Product Description Form:** Enzymatic Cleaner

**Application:** All fill sizes, see [www.pattersondental.com](http://www.pattersondental.com) for catalog numbers and safety data sheets.

**Company:** Patterson Dental Supply, Inc.  
1031 Mendota Heights Road  
Saint Paul, MN 55120  
800-873-7683

**Date of Disclosure**

September 21, 2021

**Ingredients**

Ingredient Name	CAS	Present on Designated List	Function (HCPA)
Water	7732-18-5	No	Solvent
Propylene Glycol	57-55-6	No	Enzyme preservative
Sodium citrate	68-04-2	No	Buffer
Protease (Subtilisin)	9014-01-1	EU Respiratory Sensitizer	Enzyme cleaning agent
Surfactant	Withheld as CBI	No	Cleaning agent
Amphoterics surfactant	Withheld as CBI	No	Cleaning agent
Blue dye	Withheld as CBI	No	Dye
Inorganic salt	Withheld as CBI	No	Extractant
Aromatic ether alcohol	Withheld as CBI		Preservative
Organic amine	Withheld as CB	No	Emulsifying agent
Diethyl phthalate	84-66-2	CA Priority Chemical CA Priority Pollutant	Fragrance Ingredient

**Contact for Additional Ingredient Information :** 800-873-7683

## **Links**

EU Respiratory Sensitizers: [https://echa.europa.eu/information-on-chemicals/cl-inventory-database?p\\_p\\_id=dissclininventory\\_WAR\\_dissclininventoryportlet&p\\_p\\_lifecycle=0&p\\_p\\_state=normal&p\\_p\\_mode=view&p\\_p\\_col\\_id=column-1&p\\_p\\_col\\_pos=1&p\\_p\\_col\\_count=2](https://echa.europa.eu/information-on-chemicals/cl-inventory-database?p_p_id=dissclininventory_WAR_dissclininventoryportlet&p_p_lifecycle=0&p_p_state=normal&p_p_mode=view&p_p_col_id=column-1&p_p_col_pos=1&p_p_col_count=2)

CA Priority Chemical

[https://biomonitoring.ca.gov/sites/default/files/downloads/PriorityChemicalsList\\_February2019.pdf](https://biomonitoring.ca.gov/sites/default/files/downloads/PriorityChemicalsList_February2019.pdf)

CA Priority Pollutant <https://www.law.cornell.edu/cfr/text/40/131.36>