






A full line of restorative solutions to fit your unique needs.

With a rich legacy of oral care products and services, 3M is one of the world's most trusted companies. 3M™ Filtek™ Dental Restoratives have been used in more than 1 billion restorations worldwide. Our family of Filtek dental restoratives lets you choose just the right option to fit any patient situation—from highly esthetic, to highly simple for an easier day.

Versatility & Esthetics		Esthetics & Simplicity	Efficiency & Simplicity	
<p>3M™ Filtek™ Supreme Ultra Universal Restorative</p> <ul style="list-style-type: none"> • Exceptional esthetics • Natural-looking fluorescence • Excellent polish retention • Exceptional handling • Outstanding strength for anterior and posterior • Wear resistance similar to enamel¹ • Can can be safely warmed² and injection molded • TRUE nanotechnology • 36 shades, 4 opacities 	<p>3M™ Filtek™ Supreme Flowable Restorative</p> <ul style="list-style-type: none"> • Virtually no bubbles or material “run-on” • Easy-to-bend cannula for improved access • Ergonomic design is easy to hold and inject • Excellent adaptation, polish retention and wear resistance • Can can be safely warmed^{2,3} and injection molded • Finer tip (27 gauge) enables precise placement in hard-to-reach areas. • 12 shades 	<p>3M™ Filtek™ Universal Restorative</p> <ul style="list-style-type: none"> • Predictable, repeatable esthetics • Long-lasting polish • Excellent, creamy handling • Low shrinkage stress/stress relief • Excellent wear resistance • Can can be safely warmed² and injection molded • TRUE nanotechnology • 8 designer shades plus an Extra White (XW) cover the 19 VITA classical and bleach shades⁴ • Pink Opaquer with excellent masking ability 	<p>3M™ Filtek™ One Bulk Fill Restorative</p> <ul style="list-style-type: none"> • Increased opacity for improved esthetics • One-step placement up to 5 mm • Excellent handling • Excellent adaptation • Stress relief • Excellent wear resistance • Can can be safely warmed² and injection molded • TRUE nanotechnology • 5 shades 	<p>3M™ Filtek™ Bulk Fill Flowable Restorative</p> <ul style="list-style-type: none"> • Virtually no bubbles or material “run-on” • Easy-to-bend cannula for improved access • Ergonomic design is easy to hold and inject • Injectable, one-step placement up to 4 mm • Excellent adaptation, polish retention and low shrinkage • Can can be safely warmed^{2,3}, and injection molded • 4 shades 

¹ In vitro data on file.

² Warm capsules to up to 70°C (158°F) for up to 1 hour.

³ Warm syringes to up to 70°C (158°F) for up to 1 hour, up to 25 times.

⁴ 3M internal data.

We have you covered for every class of direct restoration.

	Versatility & Esthetics		Esthetics & Simplicity	Efficiency & Simplicity	
	3M™ Filtek™ Supreme Ultra Universal Restorative (Incremental up to 2 mm)	3M™ Filtek™ Supreme Flowable Restorative (Incremental up to 2 mm)	3M™ Filtek™ Universal Restorative (Incremental up to 2 mm)	3M™ Filtek™ One Bulk Fill Restorative (Bulk up to 5 mm)	3M™ Filtek™ Bulk Fill Flowable Restorative (Incremental up to 4 mm)
Class I 	✓	✓	✓	✓ (up to 4 mm)	✓ (Liner/Base)
Class II 	✓	✓	✓	✓ (up to 5 mm)	✓ (Liner/Base)
Class III 	✓	✓	✓		✓
Class IV 	✓	✓	✓		
Class V 	✓	✓	✓		✓
Minimally Invasive	✓	✓	✓	✓	✓
Stress Bearing	✓	✓	✓	✓	

✓ = recommended use

Refer to Instructions for Use for full list of indications.



www.3M.com/Filtek

3M Oral Care
 2510 Conway Avenue
 St. Paul, MN 55144-1000 USA
 Phone 1-800-634-2249
 Web 3M.com/dental

3M and Filtek are trademarks of 3M.
 © 3M 2022. All rights reserved.
 (REV. 10/2022)
 70-2013-0807-2