


## Midwest<sup>®</sup> 45

Maximum access.  
Maximum visibility.



## The high-speed air handpiece that reaches where others can't.

### Key Benefits of the new Midwest 45

- Designed for procedures in hard-to-reach areas, such as sectioning and removal of third molars
- Compact 45° head for maximum access and visibility
- Safe, rear-facing exhaust vents direct air flow away from the surgical site, reducing the risk of air embolism
- Excellent contamination control with zero back suction
- High torque — up to 20 watts of power

### Midwest<sup>®</sup> 45

Bur Compatibility	Friction Grip 25mm length
Connection types	Fixed Backend and K-Style Coupler
Water spray	Three ports
Chuck style	Push button
Head size	10.5mm W x 13.2mm H
Weight	<60 grams
Warranty	1 year

475044  
Midwest 45  
Fixed Backend



475244  
Midwest 45  
Coupler End





## Midwest® Carbide & Diamond Burs

### The #1 Choice of Dentists

Midwest® Carbide and Diamond Burs are America's leading brand, providing the highest quality rotary cutting instruments to the dental industry for half a century, with an extensive selection of shapes, diameters and shank configurations.

To easily manage your bur inventory, download the **Midwest Bur Advisor**, a virtual bur wallet and utilities program. It's free at the Apple Store, on Google Play, or through [www.MidwestBurApp.com](http://www.MidwestBurApp.com).

### The automated maintenance system for optimal performance.

In clinical studies, correct cleaning and lubrication have proven to be the most important factors for handpiece performance and longevity. To ensure you get the most out of your handpieces, the Midwest Automate system:

- Reduces handling time by more than 80%

## Midwest® Automate®

- Automatically removes debris, precisely lubricates all moving parts, and accurately expels excess
- Ensures consistent maintenance every time
- Ends the excess mess of manual handpiece maintenance
- Optimizes long-term performance with quality maintenance



## Order Information

### Midwest® 45 High-Speed Air Handpieces

*Non-fiber optic*

*1-year warranty protection*

745044 Midwest® 45 Fixed Backend

745244 Midwest® 45 Coupler End

### Midwest® Automate® Automated Maintenance Unit

380200 Midwest Automate Unit

380520 Absorption Pad (6 Pack)

380180 Midwest Plus™ Automate Spray



Midwest Customer Service  
1-800-800-2888  
[www.dentsplysirona.com](http://www.dentsplysirona.com)