



**Product Name:** Patterson® pdCARE™ Sterilizing and Disinfecting Solution, 3.4% Glutaraldehyde, 1 Gallon Bottle

**Product Description Form:** High Level Disinfectant

**Application:** All fill sizes, see [www.pattersondental.com](http://www.pattersondental.com) for catalog numbers and safety data sheets.

**Company:** Patterson Dental Supply, Inc.  
1031 Mendota Heights Road  
Saint Paul, MN 55120  
800-873-7683

**Date of Disclosure**

December 8, 2021

**Ingredients**

Ingredient Name	CAS	Present on Designated List	Function (HCPA)
Water	7732-18-5	No	Solvent
Glutaraldehyde	111-30-8	EU Respiratory Sensitizer OEHHA REL ; CA Non-Cancer Hazards	Active ingredient
Surfactant	Withheld as CBI	No	Cleaning agent
Emulsifier	Withheld as CBI	No	Emulsifier
Anionic Surfactants	Withheld as CBI	No	Cleaning Agent
Peppermint Oil	68917-18-0	No	Fragrance Ingredient
Acetic Acid	64-19-7	No	pH Adjuster or Stabilizer

**Contact for Additional Ingredient Information:** 800-873-7683

**Links**

EU Respiratory Sensitizers: [https://echa.europa.eu/information-on-chemicals/cl-inventory-database?p\\_p\\_id=dissclininventory\\_WAR\\_dissclininventoryportlet&p\\_p\\_lifecycle=0&p\\_p\\_state=normal&p\\_p\\_mode=view&p\\_p\\_col\\_id=column-1&p\\_p\\_col\\_pos=1&p\\_p\\_col\\_count=2](https://echa.europa.eu/information-on-chemicals/cl-inventory-database?p_p_id=dissclininventory_WAR_dissclininventoryportlet&p_p_lifecycle=0&p_p_state=normal&p_p_mode=view&p_p_col_id=column-1&p_p_col_pos=1&p_p_col_count=2)

OEHHA REL ; CA Non-Cancer Hazards: <https://oehha.ca.gov/air/general-info/oehha-acute-8-hour-and-chronic-reference-exposure-level-rel-summary>