

Maxso® EP12-NS

Hygiene Handpieces

Beyes EP12-NS hygiene handpiece is a 4:1 reduction straight nose cone designed to be used with disposable prophylaxis angle for polishing. The EP12-NS offers the highest quality ceramic bearings for long-lasting smooth operation and a slim, lightweight, and compact handpiece design for ultimate user comfort.

- Beyes EP12-NS will fit all standard ISO E-Type air motor or electric micromotor
- 4:1 reduction straight nose
- Ceramic bearings
- Amazingly quiet operation
- Autoclavable
- 1 Year warranty

Order Info

Code: SL1002
Model: EP12-NS



For Disposable
Prophy Angles



Key Specification

Transmission ratio	4:1
Connection	ISO E-Type
For use with	Disposable prophylaxis angles
Speed type	Low speed
Maximum speed	10,000 rpm
Optic	Non optic
Overall length	64 mm
Shape	Straight
Weight	57 g
Sterilization	Autoclavable up to 135°C
Warranty	1 year warranty

Beyes®

The Next Standard

Beyes Dental Canada Inc.
23-595 Middlefield Road
Toronto, Ontario, M1V 3S2
Canada
Tel: 1-416-281-1888
Email: info@beyes.ca
Website: www.beyes.ca

Your Dealer for Beyes Products

Scan the QR code for more information on EP12-NS



Printed in Canada
Document # ENB021
Rev.2/21.08.19

© Beyes Canada Inc. 2019. All Rights Reserved.
Specifications are subject to change without notice.
All brands are holders of their respective trademarks.
All Beyes products are manufactured in strict compliance with Health Canada regulations ensuring quality.
Manufacturing facility is ISO 13485 certified under Canadian Medical Devices Conformity Assessment System (CMDCAS).