



IMN41637CF
16 Instrument IMS™ Infinity Series™ Cassette, Black and Orange

IMS™ Cassettes are made of long-lasting, low-maintenance stainless steel. Soft, colored silicone rails securely hold instruments throughout the processing cycle. Available in over 100 different sizes, configurations and color options, and can be specialized to meet your needs.

Why You Should Use an Instrument Management System (IMS):

By standardizing your instruments with IMS, additional time and revenue is generated while extending the useful life of your Hu-Friedy instruments. Dental professionals state that by implementing IMS, they save 5 to 10 minutes per procedure allowing the staff to spend more time on value-added activities*. In addition, safer transportation of soiled instruments reduces risk of injury and exposure to blood borne pathogens.

* When compared to single instrument decontamination. Based on market survey results.



DISCOVER
 MORE AT IMS
 WEBPAGE



For more information visit [HuFriedyGroup.com](https://www.hufriedygroup.com)

Hu-Friedy Mfg. Co., LLC, 1666 E. Touhy Ave., Des Plaines, IL 60018 | [HuFriedyGroup.com](https://www.hufriedygroup.com)
 All company and product names are trademarks of Hu-Friedy Mfg. Co., LLC, its affiliates or related companies, unless otherwise noted.
 Marks not registered in all jurisdictions.
 ©2022 Hu-Friedy Mfg. Co., LLC. All rights reserved. HFL-761/1122



HuFriedyGroup Key Opinion Leader

Dr. Shayan Barootchi & Dr. Lorenzo Tavelli

Tunneled Coronally Advanced Flap Kit



DR. SHAYAN BAROOTCHI

- Focus on: **Periodontal and Peri-implant Regenerative Surgery**
- Country: **USA**
- Favorite Instrument: **NHDPV**
- University: **University of Michigan**
- Inspired by: **Prof. William Giannobile**

DR. LORENZO TAVELLI

- Focus on: **Periodontal and Peri-implant Soft Tissue Reconstruction**
- Country: **USA**
- Favorite Instrument: **NHDPV**
- University: **Harvard School of Dental Medicine**
- Inspired by: **Prof. Giulio Rasperini**

Grafting With Science

“
 Innovation comes with persistence and curiosity.

Dr. Shayan Barootchi ”

“
 It is just a detail, but it makes all the difference.

Dr. Lorenzo Tavelli ”


TUNNELED CORONALLY ADVANCED FLAP KIT



TUNNELED CORONALLY ADVANCED FLAP KIT

TCAFKIT

Tunneled Coronally Advanced Flap Kit



Exclusive Laser Marking
for Quick Identification



TKN1

Tunneling Knife 1

For incisions that remove or recontour soft tissue. Also useful to excise interproximal tissue. Featuring a tunnelling knife.



TKN2

Tunneling Knife 2

For incisions that remove or recontour soft tissue. Also useful to excise interproximal tissue. Featuring a tunnelling knife.



PMNANOSS

Pink Matters® Tunneling Knife

Designed to help in the creation of a partial or full thickness envelope flap during connective tissue grafting procedures in perio-plastic surgery and implant surgery.



PPBUH2MX

Modified Buser/Hourigan Periosteal, Large, Black Line

Our innovative Black Line Periosteal is designed with a matte handle and performance engineered coating to optimize performance and precision.



SP20SLEM

Micro Corn Suture Pliers, Diamond Dusted

18cm. Ideal for use with guided tissue regeneration procedures.



SHDPV

Microsurgical Scalpel Handle, Stainless Steel

Mini Blade Scalpel handle used for cutting gingival tissue and making surgical incisions.



10-130-05E

5 European Scalpel Handle

This Round Scalpel Blade Handle is used for cutting gingival tissue and making surgical incisions.



DPU17

17 Dressing Pliers (Utility Pick-Up)

Made with stainless steel & feature a serrated tip for a positive grip & locking mechanism..



MIR5HDB

5 Single-Sided HD Black Line Mirror, Cone Socket

Offers superior visibility throughout any dental procedure. The Diamond Like Carbon (DLC) coating of the handle and mirror frame reduces glare up to 80%.



MH6B

Cone Socket Black Line Mirror Handle, Satin Steel™ Alloy

Offers superior visibility throughout any dental procedure. The Diamond Like Carbon (DLC) coating of the handle and mirror frame reduces glare up to 80%.



PCPUNC15CH

Qulix UNC Probe, 1-15, SE

Single-Ended, Color-Coded Probe. Ergonomic handle. University of North Carolina.



SAS1/29E2

Gracey 1/2 Mini Five™ Gracey Curette, #9 Handle

Longer terminal shank, standard diameter. Double-ended handle.



SPSPV

Microsurgical Scissor Steel, Swiss Perio

Microsurgical Castroviejo Scissors with sharp, slightly curved working ends made for fine resections of soft tissues and trimming wound margins or mucosal grafts.



NHDPV

Microsurgical Needle Holder

Microsurgical Needle Holders are used to grasp and guide the needle when suturing.



TPDAPV

Microsurgical Anatomical Tissue Pliers, Straight, Steel


Microsurgical Tissue Pliers used to grasp and stabilize soft tissue, 18mm



S18SC

Iris Super Cut Curved Scissors

Iris Curved Super-Cut scissors with 1 serrated blade used for trimming tissue or cutting sutures.



ALLENCARD

Allen Membrane Measurement Card

The Allen Membrane Measurement Card is used as a study, sterilizable surface to measure and trim membrane material.