

Fashioned for fit and performance.



Introducing  
**bluephase<sup>®</sup>**  
*style*

LED CURING LIGHT


  
ivoclar  
vivadent<sup>®</sup>  
passion vision innovation

# Slim. Powerful. Precise.

Introducing the latest innovation in LED curing light technology in a slim and lightweight design.

## Polywave™ Technology

bluephase® style's Polywave™ technology provides a broadband spectrum of 385-515 nm, which allows you to cure all dental materials on the market today!<sup>1</sup>



## Intuitive Operation

The easy-to-use 2-button operation makes the bluephase® style extremely user friendly!



## Click & Cure!

The innovative Click & Cure corded battery backup option enables you to go from cordless to corded operation in just one "click"!



## Ordering Information:

635153BU    bluephase style 100-240V  
1 – handpiece  
1 – 10mm light probe  
1 – battery  
1 – charging base  
1 – power supply  
50 – protective sleeves

## Newly Designed Light Probe

The newly designed light probe maximizes posterior access, while the 10 mm diameter of the probe allows you to cure even the largest restorations in just 10 seconds.



## Continuous Curing

Due to an energy-efficient LED technology, bluephase® style produces less internal heat compared to other LED lights, providing continuous curing with NO clinical limitations!



## Inductive Charging

The inductive charging system eliminates the need for battery contacts, which allows for much easier maintenance and charging of your battery!



**bluephase®**  
**style**  
LED CURING LIGHT

100% CUSTOMER SATISFACTION  
GUARANTEED!

[ivoclarvivadent.com](http://ivoclarvivadent.com)

Call us toll free at 1-800-533-6825 in the U.S., 1-800-263-8182 in Canada.  
©2012 Ivoclar Vivadent, Inc. bluephase and polywave are trademarks of Ivoclar Vivadent.  
<sup>1</sup>All materials known to date in the range of 385 – 515 nm  
642544 Rev 1 2/2012 4521

**ivoclar**  
**vivadent**  
passion vision innovation