

MIDWEST®

SKILLFUL HANDLING. POWERFUL DRIVE.

Introducing MIDWEST® *E* – The electric handpiece system that provides a necessary balance of power and control for all your high-speed, low-speed, and endodontic procedures

NEW

MIDWEST® *E*
Electric Handpiece System

NEW MIDWEST® E
Electric Handpiece System

THE NECESSARY BALANCE OF POWER AND CONTROL

E Power—lightweight autoclavable motor weighs only 59 g—among the lightest offered

E Control—touchscreen user interface, small handpiece designed for easy maneuvering, and the ability to make a variety of adjustments on a single screen

E Attachments offered in high speed, low speed and endodontic configurations

E Backed by the Legendary Midwest quality and support



ULTIMATE CONTROL

- Touch screen control offers one-screen review of operating settings
- Endo functionality includes user defined torque limits and 3 auto-reverse options

ATTACHMENTS

SYSTEM COMPONENTS

SKU #	ITEM
876000	Midwest® Electric Handpiece System Includes: CONTROL UNIT, MOTOR, POWER SUPPLY AND HOSE
875005	Midwest® E High Speed Attachment 1:5 HS FO CA (1:5L)
875010	Midwest® E Attachment 1:1 LS FO CA (1:1L)
875015	Midwest® E Attachment 1:1 LS CA (1:1)
875020	Midwest® E Speed Reducing Attachment 5:1 LS FO CA (5:1L)
875025	Midwest® E Speed Reducing Attachment 5:1 LS CA (5:1)
875030	Midwest® E Straight Attachment 1:1 Straight (1:1S)
875040	Midwest® E Micro Electric Control Unit Only
875045	Midwest® E Micro Electric Motor Only
875050	Midwest® E Micro Electric Motor Hose Only
875080	Midwest® E Remote Mounting Kit

For more information or to request a free demo, call **1.800.989.8825** or visit www.dentsply.com.

