



## IMPLANT 20:1 Powerful to the smallest detail.

- [HIGHLIGHTS](#)
- [PROPERTIES](#)
- [TECHNICAL DATA](#)

### Technical Data

Weight	72 g
Reduction ratio	20:1
Max. working speed	2,000 rpm*
Max. torque	70 Ncm
Head height w/o bur	13.8 mm
Head diameter	10.3 mm
Fittings	Ø 2.35 mm CA with internal cooling
Chuck	Push button
Motor interface	ISO
REF	60 43 702

\* In combination with an electric motor or with a motor speed of 40,000 rpm.