



T2 Comfort-Class Quality you can rely on.

- [HIGHLIGHTS](#)
- [PROPERTIES](#)
- [TECHNICAL DATA](#)
- [DOWNLOADS](#)

Technical data

Technical data	T2 CONTROL	T2 Boost	T2 mini
Power (W)	Up to 22	Up to 23	Up to 20
Non-load speed (rpm)	Approx. 250.000	Approx. 370.000	Approx. 400.000
Head diameter (mm)	Approx. 12,4	Approx. 11,9	Approx. 10,2
Head height (mm)	Approx. 13,8	Approx. 13,0	Approx. 12,7
Light/illuminance (lux)	Approx. 25.000	Approx. 25.000	Approx. 25.000
Wt. of high-speed handpiece (g)	Approx. 57	Approx. 56	Approx. 55
Recommended burr size (ISO)	314/315	314/315	313
Driving air pressure (bar)	3,0	3,0	3,0
Warranty (years)	2	2	2
Interface type	T2 CONTROL	T2 Boost	T2 mini
	REF	REF	REF
Sirona	63 86 648	63 86 655	63 86 663
KaVo	63 87 737	63 87 778	63 87 802
W&H	63 87 752	63 87 786	63 87 810

Interface type	T2 CONTROL	T2 Boost	T2 mini
	REF	REF	REF
NSK	63 87 760	63 87 794	63 87 828