



T3 Economy-Class Consistency that pays off.

- [HIGHLIGHTS](#)
- [PROPERTIES](#)
- [TECHNICAL DATA](#)
- [DOWNLOADS](#)

T3 Boost and T3 mini with Click&Go

Technical Data

	T3 Boost	T3 mini
Power (W)	Up to 23	Up to 20
Non-load speed (rpm)	approx. 370,000	approx. 400,000
Head diameter (mm)	11.9	10.2
Head height (mm)	13.0	12.7
Light/illuminance (lux)	–	–
Wt. of high-speed handpiece (g)	approx. 56	approx. 55
Recommended burr size (ISO)	314/315	313
Driving air pressure (bar)	2.7 – 3.0	2.7 – 3.0
Warranty (years)	1	1

Interface type	T3 Boost REF	T3 mini REF			
			Sirona	63 87 919	63 87 950
W&H	63 87 935	63 87 976			
NSK*	63 87 943	63 87 984			
KaVo	63 87 927	63 87 968			

T3 Racer and T3 mini with fixed connection

Technical Data

	T3 Racer	T3 mini			
Power (W)	Up to 20	Up to 17			
Non-load speed (rpm)	approx. 400,000	approx. 420,000			
Head diameter (mm)	11.4	10.2			
Head height (mm)	14.8	13.9			
Light/illuminance (lux)	–	–			
Wt. of high-speed handpiece (g)	approx. 53	approx. 52			
Recommended burr size (ISO)	314/315	313			
Driving air pressure (bar)	2.3 – 2.7	2.3 – 2.7			
Warranty (years)	1	1			
Interface type	T3 Racer REF	T3 mini REF			
			Midwest/ ISO	59 37	59 37

Interface type

**T3
Racer**

**T3
mini REF**

4 hole

268

276

Borden 2/3 hole

59 50
105

59 50 113