


BEAUTIFIL-Bulk

Flowable

Advanced Giomer Restorative





BEAUTIFIL-Bulk *Flowable*


Advanced Giomer Restorative

The newest member of the Shofu Giomer family, delivers a superior bulk fill composite with ideal opaque shade stability, low shrinkage stress and full polymerization at 4mm depth of cure. Beautifil-Bulk Flowable is the advanced Giomer restorative you can trust.

Beautifil-Bulk Flowable incorporates the unique, bioactive benefits of Giomer, a surface pre-reacted glass filler backed by 13 years of clinical success*, within a fast and easy bulk fill composite. Beautifil-Bulk Flowable is designed to overcome the challenges of current bulk fill products, with nice aesthetics, lower shrinkage stress and full polymerization.

*Gordan, V. et al, ed. IADR General Session, Seattle, Washington State Convention Center, March 23 2013. Clinical Evaluation of a Giomer Restorative System: Thirteen-Year Recall. Poster Session.

- High filler (72.5 wt%) and varied Giomer particle sizes reduce polymerization shrinkage and shrinkage stress while increasing compressive and flexural strength
- Low postoperative sensitivity
- Unique light-diffusion properties ensure complete polymerization at 4mm depth of cure
- Optimal opacity to match surrounding tooth structure
- Superior fluoride release and rechargability
- Self-leveling to increase cavity adaptation



Visit www.shofu.com or call **800.827.4638**



Order Information	
Syringe Refills, 2.4g	
Shades	Item No.
Universal	2030
Dentin	2031



Shofu Dental Corporation
1225 Stone Drive
San Marcos, CA
92078-4059, USA

SN197-0114

Advanced Giomer Restorative