

# MICROFLEX®

## ▲ SENSATION®



BLUE  
NITRILE  
NON-STERILE  
AMBIDEXTROUS  
TEXTURED FINGERS  
BEADED CUFF  
FOR SINGLE USE ONLY  
NON-LATEX

### POWDER-FREE EXAMINATION GLOVES

These extra-soft gloves provide exceptional fit and feel with advanced tactile sensitivity and comfort. Textured fingertips for an excellent grip. Ideal for use with small instruments.

- Enhanced comfort and tactile sensitivity
- Textured fingers for a great grip
- FDA cleared “Low Dermatitis” claim, which means Sensation® gloves have been clinically shown to reduce incidence of the redness and itching associated with contact dermatitis



### LOW DERMATITIS POTENTIAL

PRODUCT CODE		PACKAGING	
XS	N730	100	Gloves/Box
S	N731	10	Boxes/Case
M	N732		
L	N733		
XL	N734		

GLOVE PROFILE		ASTM D 6319	SENSATION®
Length	<i>mm (inches)</i>	min. 220 (8.7)	240 (9.5)
Cuff Thickness	<i>mm (mils)</i>	n/a	0.05 (2.0)
Palm Thickness	<i>mm (mils)</i>	min. 0.05 (2.0)	0.06 (2.4)
Finger Thickness	<i>mm (mils)</i>	min. 0.05 (2.0)	0.10 (3.9)
Tensile Strength	<i>Before Aging</i>	min. 14	29
	<i>After Aging</i>	min. 14	14
Elasticity/Elongation	<i>Before Aging</i>	min. 500	500
	<i>After Aging</i>	min. 400	400
Pinhole Defect Rate*	<i>(AQL)</i>	max. 2.5	1.5

\* Pinhole Defect Rate is FDA minimum requirement.



RECOMMENDED INDUSTRIES