



FEATURES AND BENIFITS OF SPRING SINGLE-USE FG DIAMOND BURS AND MULTI-USE TURBO BURS



Full report available upon request

- Spring Health Products Diamond Burs were Awarded" Best Products for 2014" by Clinicians Report (CR) in Their December 2014 Issue.
- Fresh Burs are Always Sharp and Fast Cutting.
- Spring Diamond Burs are The Highest Performance to Value Ratio Diamonds.
- Spring uses All Natural Diamonds that are Uniform in Size and Shape and have Excellent Diamond Coverage at The Tip Ensures Less Heat and a Faster Cut, Thus Less Trauma to The Tooth. Spring excels at covering the tip of the bur with diamonds, -the most difficult part of the bur to plate.
- Concentric Bur Shanks that are Micro-Precisely Calibrated and Use of 420 Heat Treated Stainless Steel Ensures Less Turbine Wear and Ease Resulting in a Smooth, Fast Cut and No Chatter
- Spring Turbo Diamonds can Reduce Chair Time by Eliminating the Use of an Additional Instrument. The Spiral Design Cuts a Smooth Axial Wall as a Single Instrument. Diamond-Free Cooling Zones That Allow the Cooling Water to Flow and Prevent Heat Generation by Bulk Reduction Results in Faster and Cooler Cutting With Less Trauma to the Tooth.
- Spring Diamonds have been recognized by Dental Product Shopper as Best Product for 2011 with a 4.5 rating.
- Made in the USA. Spring has been Manufacturing Diamond Burs for More Than 50 Years!

A Manchester, NH, evaluator



wrote that the:
"Spring Single-patient Diamond

felt like I was using an electric hand piece but I wasn't."

Fast Cutting, Reduces Chair Time, True Value

For a Free Sample, Spring Diamond Bur please call or email us at the number below



SPRING HEALTH PRODUCTS, INC

705 General Washington Blvd. Suite 701, Norristown, PA 19403 USA

(800)800-1680

customer@springhealthproducts.com | www.springhealthproducts.com