

MM-Seal™

Root canal sealer

MM-Seal™ is a high quality, epoxy resin-based paste/paste sealer for permanent filling of root canals. Its outstanding chemical and physical properties provide excellent sealing of root canals. Eugenol-free, biocompatible and radio-opaque. It can be used for all gutta percha obturation techniques.



- Dual-syringe design for accurate dispensing.
- Homogeneous, free of air bubbles and easy to mix.
- Penetrates superbly into the smallest lateral canals.
- Excellent biocompatibility.
- Can be worked for 35 minutes at 23°C.
- Sets in 45 minutes at 37°C, allowing stress-free placement of gutta percha cones.



MM-EDTA Cream

EDTA Cream

MM-EDTA Cream ensures effective cleaning and shaping of the root canal system. Its lubricating action smoothes file penetration and facilitates the removal of debris.



- Facilitates the cleaning and shaping of the root canal system.
- Lubricates.
- Facilitates instrument access and penetration.