

One Shape® Obturation


How to perfectly obturate radicular canals after shaping with One Shape®?

MICRO-MEGA® developed reliable and practical tools to simply obturate after shaping with One Shape®.


One Shape® GP Points



- Radio-opaque.
- Easy to introduce.
- Flexible for perfect adjustment to all root canal anatomies without bending.
- Easy-to-identify thanks to the coloured tip.
- Special taper adapted to One Shape®.



60 points per box. L29 mm.



60 points per box. L29 mm.

One Shape® Paper Points



- Quick and efficient drying.
- Shaped for easy and reliable penetration to the apex.
- Easy-to-identify thanks to the coloured tip.