

From Dirty to Clean and Everything In Between... It's BrandMax Clean!®



Ultrasonic Cleaners

Tri-Clean ultrasonic cleaners by BrandMax provide superior cleaning power which results in cleaner instruments. All Tri-Clean ultrasonic units come standard with heating element, cycle time controller and stainless steel basket, there is no need to pay extra for these valuable options when *BrandMax provides them for free!*

Take the *Tri-Clean* Challenge . . .

FEATURE	<i>Tri-Clean</i>	Competitor	BENEFIT
All stainless steel construction	✓	?	Scratch resistant, durable and easy to clean
Stainless steel drain valve	✓	?	Easy and effective draining of soiled cleaning solution
Digital timer	✓	?	Control time for worry free cleaning cycles
Superior acoustic pulse generators	✓	?	Powerful cavitations loosen bio-burden from instruments
Heating element included	✓	?	Expensive option included for free!
Stainless steel basket included	✓	?	Expensive option included for free!
Works with all popular cassette systems	✓	?	Supports instrument management
1 year warranty on unit & electronics / 2 year warranty on transducers	✓	?	Comprehensive protection for your investment
Superior customer service	✓	?	"Best in Class" product support

Model #	Cleaning Tank Capacity Gallons/Liters	Temp Range	Timer Setting	Power Sonic/Heat	Tank Dimensions L x W x H (in.)	Basket Dimensions L x W x H (in.)	Overall Dimensions L x W x H (in.)	Unit Weight (lbs.)	Shipping Weight (lbs.)	Shipping Dimensions L x W x H (in.)
U-5LH	1.05 gal / 4L	68-120 °F 20-60 °C	1 - 99 min.	120W / 200 W	9.5 x 5.5 x 6 in.	8.25 x 4.25 x 4.25 in.	11 x 7 x 10.75 in.	11 lbs.	13 lbs.	15 x 10 x 13 in.
U-5LHREC	1.05 gal / 4L	68-120 °F 20-60 °C	1 - 99 min.	120W / 200 W	9.5 x 5.5 x 6 in.	8.25 x 4.25 x 4.25 in.	16.625 x 12.625 x 11 in.	15 lbs.	25 lbs	21 x 16 x 15 in.
U-10LH	2.64 gal / 10L	68-120 °F 20-60 °C	1 - 99 min.	220W / 350 W	12 x 9.5 x 6 in.	10.5 x 8.25 x 4.25 in.	13.5 x 11 x 11.5 in.	17 lbs.	22 lbs.	18 x 15 x 15 in.
U-10LHREC	2.64 gal / 10L	68-120 °F 20-60 °C	1 - 99 min.	220W / 350 W	12 x 9.5 x 6 in.	10.5 x 8.25 x 4.25 in.	18.6 x 16.2 x 11.5 in.	26 lbs.	34 lbs.	21 x 19 x 15 in.
U-13LH	3.43 gal / 13L	68-120 °F 20-60 °C	1 - 99 min.	300W / 450 W	13 x 11.8 x 6 in.	11.6 x 10.6 x 4.25 in.	14.75 x 13.75 x 13 in.	25 lbs.	32 lbs.	20 x 17 x 15 in.
U-13LHREC	3.43 gal / 13L	68-120 °F 20-60 °C	1 - 99 min.	300W / 450W	13 x 11.8 x 6 in.	11.6 x 10.6 x 4.25 in.	19.25 x 18.10 x 13 in.	33 lbs.	41 lbs.	23 x 21 x 15 in.
U-20LH	5.06 gal / 20L	68-120 °F 20-60 °C	1 - 99 min.	430W / 650 W	20 x 12 x 6 in.	18.85 x 11.25 x 4.25 in.*	21.75 x 13.75 x 13 in.	36 lbs.	42 lbs.	28 x 17 x 15 in.
U-20LHREC	5.06 gal / 20L	68-120 °F 20-60 °C	1 - 99 min.	430W / 650 W	20 x 12 x 6 in.	18.85 x 11.25 x 4.25 in.*	26.75 x 18.6 x 13 in.	44 lbs.	53 lbs.	29 x 21 x 15 in.

* New 20L basket size, SN 3498 and higher

Supporting Oral Health For Better Health.™

Call Toll Free: 1-844-874-2532 or 404-418-8499
Email: CustomerService@BrandMax.us

www.BrandMax.us

BrandMax