

V-RAY Portable LED Headlight System

High Intensity Headlight (9,000 ft. candles)

PORTABLE

The V-Ray's (Product # H-VRAY-1) outstanding edge-to-edge clarity, pure white beam, and hands-free operation allows you to focus on your procedures, while the lightweight lithium-polymer battery provides power for all day use.

BRIGHT

With a powerful intensity of over 9,000 foot candles (97,200 lux) the V-RAY headlight will be the brightest light you will ever need.

LIGHTWEIGHT

With all of our portable headlights you'll experience extremely bright, shadow-free, long-lasting illumination with unrestricted mobility that enables you to work with precision while enjoying complete comfort with a headlight weight under 0.7 ounces (20g).




**9,000 ft./candles of intensity
from under 0.7 oz / (20 g)**

COMPARE THE BEAM OF THE V-RAY


Note the V-RAY's edge-to-edge consistency and uniformity



V-RAY



Fiber Optic



Other LED

ULTRALIGHT
EXTREME POWER
LONG LIFE

Lithium Polymer
Battery Technology

Mounting Options



Loupe Mount
(For All Major
Manufacturers)



Surgical
Headband



UltraLight
Headband



Del Rey Safety
Frame (Loupe-Free)

The V-Ray System (H-VRAY-1) Includes:

- Ultra high-intensity LED light source, weighing less than one ounce. Virtually maintenance free.
- One UltraLight Lithium Polymer Gel 6+ hour battery with charger (120-240V). ETL/CE UN T1-T8 Approved.
- Mounting clip available for SheerVision and most leading brands of loupes, or choice of headband or frame.
- A loupe-free frame mount is also available as shown above, with a Rx lens option for the Del Rey.
- Storage case, user manual, and a one year warranty.
- Composite FlipFilter to prevent the curing of dental composite materials.



FlipFilter



SheerVision, Inc.
4030 Palos Verdes Drive North
Suite 104
Rolling Hills Estates, California 90274
United States of America

US/Canada Toll Free: (877) 678-4274
International: + 1 310 265 8918
Fax: (310) 265-8919
customerservice@sheervision.com
www.sheervision.com

As an industry innovator, SheerVision[®] has earned a reputation for design leadership and value in optical technology supporting dental, surgical, and veterinary professionals worldwide.

SheerVision is a public company, with our stock trading under the SVSO.PK symbol.

THE LIGHTEST HIGH-PERFORMANCE LOUPES IN THE GAME