



WWW.TIPPIMMAGING.COM

Designed to FIT most digital sensors

See sizing chart below

Small Style FDS #100 Dimensions: 7.465"L x 1.638"W Head size: 1.262"

BRAND	KODAK 6100	SCHICK ELITE	SCHICK	GENDEX GX-S	EVA	XDR	INSTRUMENTARIUM SORDEX	PROGENY	MYRAY	PLANMECA	DR SUNI PLUS
SIZE	#0	#0	#0 & #1	#1	#1	#1	#1	#2	#1	#0	#0

Medium Style FDS #200 Dimensions: 7.465"L x 1.638"W Head size: 1.342"

BRAND	KODAK 6100	SCHICK ELITE	SCHICK	GENDEX GX-S	EVA	XDR	PLANMECA	DR SUNI PLUS	SUNIRAY	SOPIX 2
SIZE	#1	#1	#2	#2	#2	#1	#1	#1	#1	#1

Large Style FDS #300 Dimensions: 7.465"L x 1.88"W Head size: 1.5065"

BRAND	KODAK 6100	SCHICK ELITE	SCHICK #33	DEXIS	XDR	EVA +	INSTRUMENTARIUM SORDEX	MYRAY	PLANMECA	DR SUNI PLUS	SUNIRAY	SOPIX 2
SIZE	#2	#2	#33		#2	#2	#2	#2	#1	#2	#2	#2