



**Works Better.  
Lasts Longer.  
Costs Less.**



# STERILIZATION



## “Works Better”

Constructed with a built-in, lead-free Internal Indicator, giving you peace of mind when it comes to the sterilization of your instruments.

## “Lasts Longer”

Eliminates the need for an internal indicator strip inside the pouch with your instruments. In addition, the internal indicator will still change color when the optimal temperature is reached within the pouch.

## “Costs Less”

Compare with other dual indicator pouches to see your savings add up. Also available in bulk sizes.

# DEFEND+PLUS Pouches – with Built-in Internal Indicator



DEFEND+PLUS pouches are constructed with a built-in, lead-free Internal Indicator, giving you peace of mind when it comes to the sterilization of your instruments. The Internal Indicator ensures the correct sterilization temperature is reached where you really need it - inside the instrument compartment. This eliminates the need for an internal indicator strip inside the pouch with your instruments. In addition, the external indicator will still change color when the optimal temperature is reached. 60g medical grade ArjoWiggins paper to prevent tears, tinted plastic to help reveal any punctures, thumb notch to assist in opening, triple seal with additional chevron for extra strength.

Packaging:	Item No:	Compares to*:
<b>200 per box</b>	<b>SP-0150 - 2.25" X 4"</b> <b>SP-0250 - 2.25" X 2.75"</b> <b>SP-0350 - 2" X 7.75"</b> <b>SP-0550 - 2.75" X 10"</b> <b>SP-1500 - 3.5" X 5.25"</b> <b>SP-2500 - 3.5" X 10"</b> <b>SP-3500 - 5.25" X 10"</b> <b>SP-4500 - 7.5" X 13"</b> <b>SP-5500 - 9" X 15"</b> <b>SP-8000 - 13" X 20"</b> <b>SP-2150 - 3.5" X 10"</b>	<b>Self-Seal Sterilization Pouches with Dual Indicator - Henry Schein® Duo-Check™ Sterilization Pouches</b>
<b>1000 per box</b>		

\*Legal disclaimers: Self-Sealing Sterilization Pouches with Dual Indicator is a registered trademark of Henry Schein®. Duo-Check™ is a registered trademark of Crosstex. Use of these trademarks does not imply any affiliation with or endorsement by these companies. Mydent International is not responsible for any inaccuracy or typographical errors in this document.

**Request FREE Samples at: [www.defend.com/freesamples](http://www.defend.com/freesamples)**



80 Suffolk Court, Hauppauge, NY 11788  
 800-275-0020 • FAX 631-434-7750  
[www.defend.com](http://www.defend.com) • EMAIL: [sales@defend.com](mailto:sales@defend.com)



Follow us on Facebook.com/MydentDefend  
 Follow us on Twitter.com/Mydent\_Defend

