

# VALUMAX<sup>®</sup> Healthcare

**We Make It Better For Less**

To order: Call your local dealer • For Product Support: Tel 484-223-3103 or 212-247-7276 • email: info@valumaxhealthcare.com

## 3-in-1 Sensitive

# ArchAway

## Double-Seal

### Ear-Loop Masks

#### Features:


- Arch-Away design - Plastic wire across middle section of masks - forms an arch to keep masks away from nose and mouth at all times
  - Most stretchable ear-loop on the market, eliminate pressure to ears
  - 3-in-1 Features:  
High filtration: PFE 99.7% at 0.1micron, BFE 99.5%  
Sensitive: Silk soft inner layer eliminate any irritation to skin;  
Breathable: Delta P <2.2
  - Highly Fluid Resistant, meets ASTM Level 2 Standard, 120 mmHg
  - Manufactured in controlled cleanroom environment
  - Low particle shedding, rated category 1 for Helmky Drum Testing
- 16 attractive colors to brighten up any operations, latex free


50 pcs/box, 12 boxes/case



 Sky Blue 5630ADE-SB	 MultiColor 5630ADE-PGO	 Lt. Pink 5630ADE-LP	 Orchid 5630ADE-OC
 AQUA 5630ADE-AQ	 Plum 5630ADE-PU	 Green 5630ADE-GR	 Fern 5630ADE-FN
 White 5630ADE-WH	 Purple 5630ADE-PP	 Teal 5630ADE-TE	 Plumeria 5630ADE-PR
 Brush 5630ADE-BR	 Raspberry 5630ADE-RB	 Yellow 5630ADE-YE	 Daisy 5630ADE-DS

 mint scented

 AQUA  
5630MADE-AQ

 Brush  
5630MADE-BR

 White  
5630MADE-WH

 Sky Blue  
5630MADE-SB

Distributed by: