



Whip Mix offers you a full selection of zirconia discs. From High-Translucency to High-Translucency Xtreme, and from unshaded to individual and group shades.



Zirconia

Whip Mix offers a wide variety of zirconia in high translucency natural white and pre-shaded discs (in Vita Classic shades A1, A2, A3, BL, B1 and B2) for anterior and posterior placement. Use Vericore® ZR HT and ZR HTX for your substructures and monolithic, high-translucency, full-contour needs.

The outstanding mechanical characteristics, excellent chemical durability and the unbeaten biocompatibility combined with natural translucency make Vericore ZR HT and ZR HTX the ideal material for use in dental milling systems.

ZR HT is the ideal material to use in a posterior application due to its extremely high strength and natural appearance. It may be made as a substructure and veneered with porcelain or designed and milled as a full-contour, monolithic crown or bridge.

ZR HTX is perfect for a full contour anterior crown or short bridge. It mirrors the translucency of lithium disilicate and combines it with high flexural strength. Create a beautiful monolithic restoration and make the patient very happy.

Whip Mix Vericore HT and HTX Zirconia discs are available pre-shaded in two ways (as indicated below). The first is the conventional VITA shades A1, A2, A3, B1, B2, BL (Bleach). They are also available in convenient, inventory-friendly Group Shades. Though not exact shade matches, these group shades are a close match to each of the shades within the group. Slight stain modification (pre- or post-sintering) brings these discs to the shade match you are trying to achieve. The grouping of these shades makes keeping a large variety of shades and sizes more feasible for some laboratories.

All discs are 98mm in diameter with collar, and are available in the following thicknesses:

		10mm	12mm	14mm	16mm	20mm	25mm
VERICORE ZR HT	 Zirconia	72801	72802	72803	72804	72806	72807
	Shaded Zirconia A1	72808	72809	72810	72811	72813	72814
	Shaded Zirconia A2	72815	72816	72817	72818	72820	72821
	Shaded Zirconia A3	72822	72823	72824	72825	72826	72827
	Shaded Zirconia BL	72842	72843	72844	72845	72847	72848
	Shaded Zirconia B1	72849	72850	72851	72852	72854	72855
	Shaded Zirconia B2	72856	72857	72858	72859	72861	72862
VERICORE ZR HT Group Shades	WM1 (A1,B1,C1)	72928	72929	72930	72931	72933	72934
	WM2 (A2,A3,B2)	72935	72936	72937	72938	72940	72941
	WM3 (A3.5,B3,C2,D2,D3,D4)	72942	72943	72944	72945	72947	72948
	WM4 (A4,B4,C3,C4)	72949	72950	72951	72952	72954	72955
VERICORE ZR HTX	 Zirconia	74101	74102	74103	74104	74106	74107
	Shaded – A1	74108	74109	74110	74111	74113	74114
	Shaded – A2	74115	74116	74117	74118	74120	74121
	Shaded – A3	74122	74123	74124	74125	74126	74127
	Shaded – BL	74142	74143	74144	74145	74147	74148
	Shaded – B1	74149	74150	74151	74152	74154	74155
	Shaded – B2	74156	74157	74158	74159	74161	74162
VERICORE ZR HTX Group Shades	WMX1 (A1, B1)	74228	74229	74230	74231	74233	74234
	WMX2 (A2,B2,C1,C2)	74235	74236	74237	74238	74240	74241
	WMX3 (A3,D2,D3)	74242	74243	74244	74245	74247	74248
	WMX4 (A 3.5, B3, B4)	74249	74250	74251	74252	74254	74255
	WMX5 (A4,C3,C4,D4)	74256	74257	74258	74259	74261	74262