

# **hi-tray light edentulous**

## *Anatomic and resistant.*

Anatomic impression trays in **green rigid plastic** for edentulous arches.



HI-TRAY LIGHT EDENTULOUS

### FEATURES

- Made in **green** plastic
- Available in the following **versions**:
  - small** (upper and lower)
  - medium** (upper and lower)
  - large** (upper and lower)

### ADVANTAGES

- **Rigid** and **non-deformable**
- **Disposable**
- **Adaptable**
- **Perforated design** for **improved mechanical retention**
- **Reinforced grip** with ribbing for **improving the distribution of the pressure** exerted



**hi-tray**  
**Aimed rigidity.**