

## MICRO-THIN FOR ADVANCED COMFORT, FIT AND FEEL

**MICRO-TOUCH Micro-Thin Nitrile** examination gloves feature an advanced stretch formulation providing ultimate comfort and fit. With its high strength and excellent protection features, this glove is ideal for use in nearly all clinical settings. Featured in a 300 count box to help you go GREEN while driving operational efficiencies. More gloves in a box means less packaging waste into the environment.

### Safety & Protection

**Infection Control**

**Allergy Concerns.** This examination glove is made of 100% nitrile. It does not contain natural rubber latex, thus addressing Type I (latex) allergy concerns.

**Sterilization.** This is a non-sterile examination glove.

**Quality Benchmarks.** Meets or exceeds ASTM examination glove standards. Manufactured within the quality guidelines of ISO 13485:2003, ISO 9001:2008 and FDA-QSR.

**Cuff.** Inverted, 9.5" beaded cuff.

**Resistance.** Nitrile gloves are recommended for conditions where high strength and chemical protection are required.

### Skin Protection

**Powder-Free.** To help eliminate powder-related complications and sensitivity.

### Comfort & Fit

**Hand Fatigue**

**Stretchability (Modulus).** Nitrile has high elasticity providing wearers with good comfort and fit.

**Shape of Glove/Former.** Ambidextrous.

### Tactile Sensitivity

**Optimal Feel.** The advanced nitrile thin-film formulation facilitates optimal 'touch and feel' to the user, instilling confidence to conduct the procedure.

**Dependability.** The textured fingertips of the glove promote sure handling of materials and instruments, in both wet and dry application environments.

**Memory.** Nitrile has high elasticity and good memory allowing the film to adapt to the wearer's hand.

### Grip

**Application Driven.** The textured fingertips of the glove provide extra protection for wet or dry grip, and promote sure handling of materials and instruments.

