

TECHNICAL DATA

MICRO-TOUCH[®] Micro-Thin Nitrile, Powder-Free, Latex-Free Medical Examination Gloves

Material 100% nitrile rubber (acrylonitrile/butadiene/methacrylic acid terpolymer).
Not made from natural rubber latex.

Physical Properties (ASTM Requirements)

	Before Aging	After Aging
Freedom from Holes (inspection level 1)	1.5 AQL*	N/A
Tensile, MPa (p.s.i.) minimum	14.0 (2,030)	14.0 (2,030)
Elongation (%) minimum	500	400

Ansell manufacturing specifications meet or exceed ASTM specifications.

Color	Blue
Accelerators and Antioxidants	Zinc Dibutyldithiocarbamate (ZDBC) Zinc Diethyldithiocarbamate (ZDEC)
Donning Agent	None (Powder-Free)
Certification	In compliance with ASTM D6319-10, D6124, D5151-06, D6978-05, D3767, D412 and F1671-07 Manufacturing facility registered ISO 9001, ISO 13485
Powder Content**	<2 mg/glove
Country of Origin	Malaysia
Packaging	300 gloves/dispenser box - 10 dispenser boxes/case
Dispenser Dimensions	245 mm x 128 mm x 83 mm
Expiry Date	Three (3) years in a cool environment, shielded from ozone, UV light, and excessive heat.
Style	603431x

Dimensional Specifications

SIZE	Average length mm (inches)	Average Palm Width mm (inches)
X-Small	230 (9.06)	80 (3.2)
Small	230 (9.06)	85 (3.4)
Medium	230 (9.06)	95 (3.8)
Large	230 (9.06)	106 (4.2)
X-Large	230 (9.06)	115 (4.5)

Thickness (Single)

Data from Sampling	Average Thickness
Fingertip, mm (mils)	0.07 (2.8)
Palm, mm (mils)	0.05 (2.0)
Cuff, mm (mils)	0.04 (1.6)

Product Information

SIZE (identified inside glove)	Reorder Numbers	ORACLE Numbers
X-Small	6034310	118255
Small	6034311	118256
Medium	6034312	118257
Large	6034313	118258
X-Large	6034314	118259



*Acceptable Quality Level - Maximum % of gloves with holes under standard sampling and testing guidelines.

**Per ASTM: maximum powder limit for powdered medical gloves is 10 mg/dm²; maximum for powder-free gloves is 2 mg/glove.

For additional information and samples, call Ansell at **800-952-9916** or visit us at www.ansellhealthcare.com

Ansell Healthcare Products LLC | ansellhealthcare.com
111 Wood Avenue South, Suite 210, Iselin, NJ, 08830, USA.

* and [™] are trademarks owned by Ansell Limited or one of its affiliates. © 2014 Ansell Limited. All rights reserved.

SS800