




TC CL II

**THE SMALLEST,
MOST POWERFUL
AND DURABLE
LED CURING
LIGHT**



Weight	2.75 oz or 83 gm
Overall Length	8" or 20.3 cm
Grip Diameter	.750" or 1.9 cm
Base Diameter	.815" or 2.1 cm
Battery/Output	45 minutes continuous run time Constant light output 1200mw/cm ²
Power	Programmable high setting 1200mw/cm ² and low setting 600mw/cm ² for optional step curing
Curing Time	5-15 seconds The CURE offers fast, reliable curing in as little as 6 seconds on most modern composites (A3 or lighter shades)
Led Wave length	440-465nm
Battery Life	Long lasting lithium-ion battery 12-24 months battery life - replaceable at chair side

Timer	Digital programmable timer - 8 settings from 2 seconds to 20 seconds
Beeper	2 second audible beeper
Housing	Durable anodized aluminum handle High-tech, light-weight all metal housing. Small in size and one of the lightest on the market. Easy to clean.
On/Off Switch	Dual switches: Ease of operation in any position.
Country of Origin	Made in the USA, including the circuit board
Warranty	1 year warranty on all parts and labor including battery



(Corded model)

THE CURE ORIGINAL LED CURING LIGHT (TC 01)

- Dual Switches: Ease of operation in any position.
- High / Low output from 600 to 1200 mw/cm².
- 16 second timed cycle, beeps every 2 seconds.
- Included with the TC-01: AC adapter power supply, 8 mm lens, lens wrench, 2 light shields, 50 protective sleeves, handpiece hanger.



Made in USA

For more information or to order, call 800.800.1680, +1.610.630.9171 or fax 610.630.9730
www.springhealthproducts.com • Product Specialist: customercare@springhealthproducts.com