

# Midwest Stylus® Plus High-Speed Handpiece

## THE NEXT GENERATION OF AIR TURBINE FROM DENTSPLY SIRONA

### More Power

All new rotor design provides up to 30% more power than previous Stylus generations.

### More Control

Dynamic speed control in the head helps maintain cutting power while also protecting against excessive speeds that can cause premature bearing wear.

### Less Noise

Unique bearing design, computer balanced rotor, and tuned air flow work together to minimize vibration and noise for a more pleasant work environment.

### Improved Design

Ceramic bearings are protected with quick-stop technology to help keep a clean service record; 2-year warranty.

792400 Midwest Stylus® Plus — Midwest backend

792410 Midwest Stylus® Plus — Sirona backend


792420 Midwest Stylus® Plus — NSK backend

792430 Midwest Stylus® Plus — K-style backend



## PERFORMANCE UNDER PRESSURE

Stay in control with the Midwest Stylus® Plus. Dynamic speed control in the head limits free speed to protect bearings against excessive RPMs at higher pressures. This also translates to lower noise without impacting peak cutting power.



## COMPLETE CUTTING SOLUTIONS

As a leader in both handpiece and bur manufacturing, only Midwest offers the ultimate integration of design, performance, and integrity.



## KEEP IT CLEAN

In clinical studies, correct cleaning and lubrication have proven to be the most important factors for handpiece performance and longevity. To ensure you get the most out of your handpieces, the Midwest® Automate® system:

- Reduces handling and cleanup time by more than 80%
- Automatically removes debris, precisely lubricates all moving parts, and accurately expels excess
- Ensures consistent maintenance every time
- Ends the excess mess of manual handpiece maintenance
- Optimizes long-term performance with quality maintenance



## ORDER INFORMATION

### Midwest Stylus® Plus High-Speed Handpieces

- 792400 Midwest Stylus® Plus — Midwest backend
- 792410 Midwest Stylus® Plus — Sirona backend
- 792420 Midwest Stylus® Plus — NSK backend
- 792430 Midwest Stylus® Plus — K-style backend

### Couplers

- 790245 5-Hole Coupler — Midwest style
- 790246 6-Pin Coupler — Midwest style
- 790335 KS5 5-Hole Coupler — K-style
- 790336 KS6 6-Pin LED Coupler — K-style

### Midwest® Automate® Automated Maintenance Unit

- 380200 Midwest® Automate® Unit
- 380090 Midwest Stylus® Maintenance Coupler
- 380450 K-style Maintenance Coupler
- 380480 NSK Maintenance Coupler
- 380510 Sirona Maintenance Coupler
- 380520 Absorption Pad (6 Pack)
- 380080 Midwest Plus™ Aerosol Spray Can



Customer Service: 1-800-800-2888  
 901 West Oakton Street  
 Des Plaines, IL 60018-1843 USA  
[www.DentsplySirona.com](http://www.DentsplySirona.com)

\*Data on file. See IFU for recommended operating pressure.