

By: TORNADO PLASTICS Inc.

TIPP ITEM#	TIPP INTRAORAL CAMERA BARRIER SOLUTIONS											
IOC100	Length 4.25"	07-110-7333										
BRAND	GENDEX	DEXIS	XDR RADIOLOGY	XDR RADIOLOGY	DÜRR DENTAL	SIRONA	SORDEX	RF AMERICA IDS	RF AMERICA IDS	RF AMERICA IDS	RF AMERICA IDS	
MODEL	C300	DEXCAM #3	XDR	Suni	VISTACAM IX	SIROCAM-3	DIGORA VIDI	DISCOVERY 360	FOCUS	DR. CAM	EZ SHOT	
IOC200	Length 3.75"	07-083-6304										
BRAND	ACTEON	ACTEON										
MODEL	SOPRO 717	SOPRO LIFE										
IOC300	Length 4.5"	07-035-0553										
BRAND	DÜRR DENTAL	AIR TECHNIQUES	AIR TECHNIQUES	AIR TECHNIQUES	DEXIS							
MODEL	VISTACAM CLIX	POLARIS	LUNA	SPECTRA	DEXCAM #3							
IOC500	Length 6.25"	07-035-0520										
BRAND	DR QUICKLOOK	SCHICK	SCHICK	SCHICK	SCHICK							
MODEL	DR QUICKLOOK	USB #1	USB #2	USB #3	USB #4							
IOC400-600	Length 7"	07-110-7358										
BRAND	MY RAY	3MESPE	KAVO	NDI	KODAK	KODAK	KODAK	SCHICK	SCHICK	SCHICK	SCHICK	
MODEL	C-U2	TRUE DEFINITION SCANNER	ERGOCAM 5	DOC PORT	CARESTREAM 1200	CARESTREAM 1500	CARESTREAM 1600	USB #1	USB #2	USB #3	USB #4	
IOC700	Length 4"	07-110-7374										
BRAND	SIRONA	SIRONA										
MODEL	SIROCAM AF	SIROCAM AF +										
IOC800	Length 8"	07-110-7382										
BRAND	AIR TECHNIQUES	AIR TECHNIQUES	AIR TECHNIQUES	AIR TECHNIQUES	DEXIS	DJRR DENTAL	DIGI-DOC IRIS #1	IRIS				
MODEL	SPECTRA G4000	Polaris	CAM X LUNA	CAM X TRITON	DEXCAM #4	VISTACAM CLIX IX	LTE	HD 3.0				
IOC900	Length 7"	07-110-7390										
BRAND	KODAK	AIR TECHNIQUES	AIR TECHNIQUES	DEXIS	ACTEON							
MODEL	CARESTREAM 1200	CAM X SPECTRA J2300	CAM X ELARA	DEXCAM #3	SOPRO 617							
NEW IOC1000	Length 5"	07-xxx-xxxx										
BRAND	DIGI-DOC											
MODEL	IRIS 2.0											
NEW IOC1100	Length 7.75"	07-xxx-xxxx										
BRAND	DIGI-DOC											
MODEL	IRIS X80											

*Certain intraoral camera brands/models may suitably fit several intraoral camera sheath styles

INTRAORAL CAMERA BARRIER SOLUTIONS (IOC) FEATURES & BENEFITS

- >>ALWAYS UPDATED FOR USE ON A VARIETY OF MOST POPULAR INTRAORAL CAMERAS
- >>SUPER CLEAR MATERIAL EXCEPTIONAL CLARITY
- >>EASILY ADJUSTABLE FITTED TIP ELIMINATES AIR POCKETS & IMAGE DISTORTION
- >>SPECIFICALLY DESIGNED FOR SEAMLESS LENS CONTACT
- >>SOFT-TOUCH MATERIAL DESIGNED FOR PATIENT COMFORT
- >>PATENTED MATERIAL COMBINATION FOR EFFORTLESS POSITIONING
- >>EASY TEAR-AWAY REMOVAL AVOIDS STICKING TO LENS & DAMAGING EQUIPMENT
- >>LATEX-FREE FOR NO UNPLEASANT FLAVOUR
- >>INTERLEAFED DISPENSER PACKAGING TO AVOIDS CROSS CONTAMINATION




TIPP BARRIER SOLUTIONS by Tornado Plastics Inc. 4805 rue Bourg, Montreal, Quebec H4T 1H9 Canada
Tel.: 1-800-745-0109 Fax: 514-748-0808 www.tippimaging.com tipp@torplast.com

CAMERA BARRIER INSERTION AND TEAR-AWAY REMOVAL INSTRUCTIONS

EASY INSERTION



TEAR-AWAY REMOVAL



- 1) While holding the bottom of sheath and paper, open and insert camera under frosted layer with lens facing paper.
- 2) Gently slide camera into sheath until tip is snugly positioned at the end of sheath.
- 3) DO NOT FORCE. If product becomes torn or punctured do not use.
- 4) Once properly positioned, detach from paper. Camera is now ready to use.
- 5) Once treatment completed, tear-away frosted layer from clear and immediately discard both.