

Patterson[®] Phosphor Plate Barriers

Instructions for Use

Safe and effective, Patterson Phosphor Plate Barriers come in sizes 0, 1, 2 and work with all dental PSP plates. All you do is load your plate, seal, expose, unload and process.

After your plate is processed and erased, re-load and you're ready to go!

Patterson Phosphor Plate Barriers do not contain Natural Rubber Latex. Intended for single use only.



LOAD
INSÉRER



SEAL
SCELLER



UNLOAD
RETIRER

Enveloppes protectrices pour plaques phosphoriques

Mode d'emploi

Sûres et efficaces, les enveloppes protectrices pour plaques phosphoriques Patterson sont offertes en tailles 0, 1 et 2, et conviennent à toutes les plaques phosphoriques dentaires. Vous n'avez qu'à insérer votre plaque, la sceller, l'exposer, la retirer et la traiter.

Une fois votre plaque traitée et effacée, insérez-la de nouveau dans une enveloppe pour recommencer le processus!

Les enveloppes protectrices pour plaques phosphoriques Patterson sont exemptes de latex de caoutchouc naturel. À usage unique seulement.

Size 0 Reorder Number
N° de produit – taille 0

129-3943

Made in China - Fabriqué en Chine

Size 1 Reorder Number
N° de produit – taille 1

129-3950

Made in China - Fabriqué en Chine

Size 2 Reorder Number
N° de produit – taille 2

129-3927

Made in China - Fabriqué en Chine