



tradeshow

Booth 2600

pocketguide

Greater New York Dental Meeting
Jacob K. Javits Convention Center
New York, NY | November 25-28, 2012



**Play
and
Win!**



› Try your luck at Patterson Dental's
trade show game

Visit the Patterson Booth to play!

The Real Deal Game

Win bonus Advantage Dollars!

Place a merchandise order at the Patterson booth and earn chips to play on The Real Deal game board. Where your chip lands determines your bonus Advantage Dollars!

Minimum payout: 100 Advantage Dollars

Maximum payout: 1,500 Advantage Dollars

Merchandise Purchases Before Tax and Shipping	Number of The Real Deal Chips Earned
\$1,000 to \$1,999	1
\$2,000 to \$2,999	2
\$3,000 to \$3,999	3
\$4,000 to \$4,999	4
\$5,000 and greater	5

Patterson Booth 2600

CAESY Cloud:

Great CAESY content is just a click away!



Connect to CAESY Cloud today!

Patterson Booth 2600

Gain instant access to more than 280 multimedia presentations featuring 3D animation, full-motion video, narration and colorful images.

With CAESY Cloud, you can access these valuable presentations from multiple locations and multiple devices throughout your office at the same time.

- New videos and features are pushed online as they become available
- Low monthly start-up and subscription fee
- No installation required – access from anywhere
- Access on multiple media devices, including iPads and iPhones
- Multiple licensing options

Patterson Booth 2600

Now add
CAESY Cloud
videos to your
website



Visit us at Booth 217!

Heraeus

NEW Venus® Temp 2 For Seamless Transitions



- Nano-particle formula for a brilliant, esthetic gloss finish
- Exceptional wear resistance and flexural strength
- Optimal flowability for easy handling and application
- Fluoride-releasing

PURCHASE 2 Venus Temp 2 Cartridges, **RECEIVE** 1 of the same **FREE!** 1000804079



Venus®
The New Esthetics

VITA Easyshade® Compact

Special GNYDM Price: 10% OFF!

The Shade Matching Technology That Reduces Remakes

Instant Shades in the Palm of Your Hand

Cordless, sleek and incredibly fast, the VITA Easyshade® Compact delivers the fastest, easiest way to measure tooth shades under any lighting condition. Easyshade Compact gives you the ability to conquer all shade challenges all day, every day, instantly and accurately. Now you can measure natural teeth, check restoration shades and track bleaching progress with the push of a button. Order today and instantly prescribe VITA 3D-Master and Classical shades.

Order now from Patterson Dental at 800-873-7683

Exclusive
Distributor

PATTERSON
DENTAL

VITA®

©Vident 2012. Easyshade and 3D-Master are registered trademarks of VITA Zahnfabrik.



Avoid bib chain contamination risk.

Use **Bib-Eze**
Disposable Bib Holders

➔ Buy 3, get 1 FREE!*



Box of 250 (226-9108) \$12.95



Patterson will ship free goods automatically with your purchase; qualifying purchase amount must be of same item number on one order.

*Offer expires: December 31, 2012.


PATTERSON
DENTAL

For more information or to place your order,
please contact your Patterson Dental
representative, local branch or call 1.800.873.7683

Trust.Worthy.Innovation.



A Breakthrough in Temporary Cement

A new generation of translucent temporary cement with color blending nanofillers

Introducing!



ZONEfree

4g, 8 mix tips (162-7801) \$44.95

LAUNCH SPECIAL!

**Buy 2 ZONEfree syringes,
get 1 FREE!**

Developed using proprietary prismatic nanofillers, ZONEfree has color adjustment potential to blend with any adjacent surface making it a perfect choice for all your temporary cementing needs.

Order ZONEfree from your Patterson dental dealer. For a list of dealers visit www.duxdental.com/dealers. Choose to be free today. Choose ZONEfree.

Patterson will ship free goods automatically with your purchase; qualifying purchase amount must be of same item number on one order.


PATTERSON
DENTAL

For more information or to place your order,
please contact your Patterson Dental
representative, local branch or call 1.800.873.7683

Photos courtesy of Shannon Pace Brinker, CDA, COD



Automix Easy Application



Effortless Clean-Up



Beautiful Esthetic Results

*Offer expires: December 31, 2012.

Trust.Worthy.Innovation.



Relief for Snoring and Sleep Apnea



Now you can offer dental sleep medicine for **your** practice with the TAP® 3 Elite. The most effective oral appliance for snoring and sleep apnea. **Period.**



- Opportunities await in the fastest growing sub-specialty of dentistry
- Incremental revenue for your practice
- Improve your patient's overall health and quality of life

TAP® Snoring and
Sleep Apnea Toolkit
Patterson Item 105-0152



Patterson Booth 2600



For more information or to place your order, contact your Patterson representative, local branch or call **1-800-873-7683**.



RelyX™ Ultimate Adhesive Resin Cement

- Industry-leading bond strength
- Fast and easy — one bottle, one syringe
- Total-etch or self-etch versatility
- Fluorescent shades for world-class esthetics
- Ideal for CAD/CAM and glass ceramic restorations

www.3MESPE.com/RelyXUltimate

Visit 3M ESPE booth #4609 to learn how ultimate bond strength can be yours with FEWER working steps!

3M, ESPE and RelyX are trademarks of 3M or 3M Deutschland GmbH.
Used under license in Canada. © 3M 2012. All rights reserved.



3M ESPE



A GREAT CHOICE

MICROBRUSH[®]
INTERNATIONAL
www.microbrush.com

VISIT US AT BOOTH 510
and receive a FREE Micro
Bear with your order!



WHEN YOU BUY



Two 4-packs of Microbrush[®]
Plus Dispenser Series or 8 tubes of
Microbrush[®] Tube Series applicators

GET YOUR CHOICE



100 Ultrabrush[®] 2.0 applicators
plus a sion[™] dispenser

OR



100 Microbrush[®] X[™] applicators

OR



50 fine size Microbrush[®] applicators
plus a sion[™] dispenser

Offer valid 10/1/12 - 12/31/12

Microbrush[®] International

866 • 866 • 8698

More Masterpieces Daily.

Everyday we are working with leading clinicians to elevate the art of dentistry. From prep to polish, we support your artistic creations with instruments created for your hands.

We are dedicated to engineering instruments that deliver efficiencies, such as **Great White**® carbides and **Piranha**® **Diamonds** to single-step **Jazz**® **Polishers**, that mean less chair time for patients and more time for you. Efficiency that is moving more and more smiles everyday.



"I have yet to experience anything that cuts as fast, accurate and smooth as a Great White® bur, saving my practice time and money while reducing stress."

- Howard S. Glazer, DDS


*Recommended choice of the Midwest
Dental Evaluation Group (MDEG)*

- Interface 2006 Newsletter



Visit us at booth #2000

Greater New York Dental Meeting | November 25 - 28 | New York



GREATER NEW YORK SHOW SPECIAL

BUY 2 KITS, GET 1 FREE*

*Of equal or lesser value



JAZZ™ Supreme
for faster, more
beautiful polishes



Great White® Ultra
a unique, faster alternative
for crown preparations



Fissurotomy®
a conservative solution
to suspect fissures



Safe End Finishing
procedure specific instruments
that deliver ultimate efficiency



Great White® Z
specifically designed
to cut through zirconia



Inlay/Onlay
for preparation of
Inlay/Onlay Restorations

To claim free goods: On or before January 15, 2013, mail a copy of your Patterson invoice to:
SS White Burs, Inc., Attn: Promotions Department, 1145 Towbin Ave., Lakewood, NJ 08701 or fax it to 732.905.0987.
Please note that these are end-user doctor promotions and must be redeemed by the qualifying practice.
Please allow 4-6 weeks for free delivery. Invoices cannot be combined for free goods redemption.

PRACTICE INSPIRATION™

SS WHITE®

800-535-2877 | www.sswhiteburs.com
1145 Towbin Avenue Lakewood, New Jersey 08701



\$129⁹⁵
 [569-3627]

FREE GIFTS

**10 Delta Silicone
 Tubes + 1 Re-Force**



Quickmat Deluxe

**The ideal sectional matrix
 system for perfect contact
 points in class II restorations**



*compare contents
 and price!*

[569-3627] Quickmat Deluxe Kit

Content: 5 Quick rings, 10 silicone tubes,
 200 contoured sectional matrices
 .0015" thin (two heights),
 100 contoured sectional matrices
 .001" microthin (two heights),
 300 assorted wood wedges.

Special Financing*

3.3%

*for 33 months on Schick
sensors, remotes, cables
and training*

*Must take delivery by December 31, 2012.

SCHICK 33

The Schick 33 sensor is the most advanced sensor on the market, delivering an unparalleled combination of high-resolution imaging and dynamic image management in a modular platform.


- Image preference presets that give users the ability to save images based on diagnostic needs, from general dentistry to endodontics, periodontics and restorative dentistry.
- User-definable tasks for multi-clinician practices make it easy to change settings according to preference or diagnostic task.
- Three sensor sizes (0, 1, 2) and three cable lengths (3', 6', 9') to choose from.



schick
by **sirona**



Preventive Products



Buy 3 Preventive Products,
Get 1 FREE!* PLUS 1 MI Paste™ Marketing Kit FREE!*
promo 1709

Choose from MI Paste™ 10-Tube Packages, MI Paste Plus™ 10-Tube Packages,
GC Dry Mouth Gel 10-Tube Packages, Saliva-Check MUTANS Kits or Saliva-Check BUFFER Kits.

Buy 2 MI Varnish™ Packages,
Get a \$25 Best Buy or Amazon Gift Card FREE!*
promo 1666

800.323.7063 • www.gcamerica.com • www.gcamerica.com/training • © 2012 GC America Inc.

*Terms and Conditions of Promotions: Effective Dates: 11/1/2012 - 11/30/2012. Promotions cannot be combined with any other promotional offers. No substitutions allowed. Offers are valid in the US and Canada only. MI Varnish & MI Paste Plus are not available in Canada. Free goods coupon must be redeemed within 6 weeks of Patterson invoice date. A copy of Patterson invoice/proof of purchase of qualifying product(s) can be faxed to 800-423-2963 or mailed to: GC America Inc., Attn: Customer Care Group, 3737 West 127th Street, Alsip, IL 60803. Free goods will be of equal or lesser value to the lowest selling unit price. Please allow 6-8 weeks for free goods delivery. All free goods will be fulfilled by GC America Inc. except where noted. Please indicate requested free goods on Patterson invoice.



Restoratives & Cements



Buy 1 GC FujiCEM™ 2 Product and 2 other Glass Ionomer Products,
Get a \$75 Best Buy or Amazon Gift Card FREE!*
 promo 1795

Excludes GC FujiCEM™ 2 Intro Kit and Dispenser

Buy 2 G-CEM Refills, **Get 1 FREE!***
 promo 1711

G-CEM™ 50-Capsule Refill:					\$296.00
205-8261 A2	208-8279 AO3	205-8295 BO1	205-8287 Translucent	205-8303 Assorted	
G-CEM™ Automix Refill:					\$134.95
205-8402 A2	205-8410 Translucent				

Visit GC America's Booth at the 2012 Greater New York Meeting
to Take Advantage of these Special Offers! (Booth 1813)



Impression Materials



Buy 3 EXA Products, **Get 1 FREE!***

promo 1708 • Excludes Putty

EXA'lence™ 4-Cartridge Regular Set Refill:				\$97.95
204-0830 Extra Light	204-0848 Light	204-0855 Monophase	204-0814 Heavy	
204-0822 Heavy Rigid				
EXA'lence™ 4-Cartridge Fast Set Refill:				\$97.95
204-0780 Extra Light	204-0798 Light	204-0806 Monophase	204-0764 Heavy	
204-0772 Heavy Rigid				
EXAFast™ 2-Cartridge Regular Set Refill:				1@ \$30.95 4@ \$29.95
326-0189 Injection	326-0205 Regular	326-0213 Monophase	326-0197 Heavy	
EXAMIX™ 2-Cartridge Fast Set Refill:				1@ \$30.95 4@ \$29.95
326-3464 Injection	326-3456 Regular	326-3480 Monophase	326-3472 Heavy	



The most clinically capable CEREC system ever.

From single crowns to multiple units, individual quadrants to full arches, CEREC Omnicam "does it all."

- Designed for Access
- Color Streaming for Confidence
- Powder Free for Convenience



Exclusive
Distributor


PATTERSON
DENTAL

sirona



BLOSSOM LATEX POWDER FREE EXAM GLOVES

DESCRIPTION	SIZE	ITEM NO.
Ambidextrous/Beaded Cuff	XS	486-4484
Non-Sterile / Disposable	S	486-4492
100 Gloves Per Box	M	486-4500
10 Boxes Per Case	L	486-4518
	XL	486-4526



BLOSSOM POWDER FREE TEXTURED EXAM GLOVES WITH ALOE VERA

DESCRIPTION	SIZE	ITEM NO.
Ambidextrous/Beaded Cuff	XS	486-3718
Non-Sterile / Disposable	S	486-3726
100 Gloves Per Box	M	486-3734
10 Boxes Per Case	L	486-3742
	XL	486-3759

SUPER
D

Ambid
Non-S
100 G
10 Bo



Caution



BLOSSOM
www.blossom-disposables.com

Supergloves Powder Free Latex Exam Gloves

DESCRIPTION	SIZE	ITEM NO.
Ambidextrous/Beaded Cuff	S	486-4112
Sterile / Disposable	M	486-4120
Gloves Per Box	L	486-4138
Boxes Per Case	XL	486-4146



Blossom Latex Powdered Exam Gloves

DESCRIPTION	SIZE	ITEM NO.
Ambidextrous/Beaded Cuff	XS	486-4344
Non-Sterile / Disposable	S	486-4310
100 Gloves Per Box	M	486-4328
10 Boxes Per Case	L	486-4336
	XL	

: "Safe use of this glove by or on latex sensitized individuals has not been established"





GROW
YOUR REVENUES

REDUCE
YOUR COSTS

WOW
YOUR PATIENTS

Revolutionize your practice with **RevenueWell**: the automated practice marketing and patient communications suite that takes the hassle out of growing your business.

- ✓ Automatically recall and reactivate patients
- ✓ Automatically confirm each appointment
- ✓ Increase your online visibility
- ✓ And so much more!

See a demo at **Patterson Dental booth 2600** or visit www.PattersonDental.com/AppStore

SHOW TIME PROMOTION!
NO Installation Fee!



pattersondental.com/AppStore



RevenueWell

**PATTERSON
DENTAL**



YOUR DATA. OUR SOFTWARE. **CHANGE MADE SIMPLE.**

Now, when you take advantage of our free practice management software,
GET A FREE DATA CONVERSION.

Transitioning to a new software system can be challenging. But Eaglesoft Practice Management Software is making it easier for you – when you switch, we'll give you your data conversion for free.* We'll walk you through it every step of the way. But hurry, this offer expires December 31, 2012.

Contact your Patterson representative or call 800.294.8504 for more information.

*Image conversion subject to standard fees.

patterson.eaglesoft.net

power *personalized*





Introducing A New "Aseptically Manufactured" Articaine

Orabloc™

(articaine HCl 4% and epinephrine 1:100,000
and epinephrine 1:200,000) injection



- Rapid onset of anesthesia within 1-9 minutes
- Complete anesthesia lasts about 1 hour for infiltration, up to 2 hours for nerve block
- Only 10% overage of epinephrine
- A 24 month shelf life at room temperature
- Very low manufacture related degradation products such as articaine acid and epinephrine sulphonic acid
- Sodium edetate free, methylparaben free and latex free
- Each Orabloc™ cartridge is sealed individually in the blister for maximum protection during storage and handling

568-2604

568-2612

Link to www.orabloc.com

For more information call Pierrel Research USA at 1-610-989-4222 or send an email to info.pharma@pierrelgroup.com



Learn more at Patterson Dental
booth 2600 and visit
www.PattersonDental.com/AppStore!



Patterson Dental Explore Store
> Apps & More!

Across multiple platforms and devices, apps and resources from Patterson's new Explore Store help you manage and grow your practice.

pattersondental.com/AppStore



PATTERSON
DENTAL



Delight your patients, delight your staff
– **enhance your practice!**

- 20% off all accessories
- \$10 off the B145 small brewer unit
- \$20 off the B150 and B155 medium-sized brewer units
- Buy one box of coffee or tea, get one box of equal or lesser value FREE (up to four boxes)

Use code VL2
when placing order.

Take advantage of our **tradeshow discounts**
– order at the Patterson booth today!

Patterson Booth 2600



PATTERSON
OFFICE SUPPLIES



PATTERSON DENTAL





©2012 Patterson Dental Supply, Inc.
All rights reserved.

P130682 (9/12)